



### Anatomy Tips

- Use the spine as a Fluoro reference for the LA. L & R borders of the spine often = L & R sides of LA.
- The Coronary Sinus (CS) is the floor of the LA.
- Normal EsoSure deflection is from the R to L spinal border in AP.

- View insertion & rotation in AP. Watching Fluoro, imagine the position of the EsoSure in 3D. Deflection occurs posterior-laterally between the lung and spine.
- Evaluate deflection from a perpendicular view: RAO for R deflection & LAO for L deflection.




### Initial R Deflection



### If deflection is poor, use the ventilator to change the anatomy.

A narrow LA to spine distance may block deflection. Fluoro shows the EsoSure curve stops at mid-spine.




### If deflection stops at mid-spine, use Valsalva

- 1) Rotate the primary curve to the desired side behind the trachea;
  - 2) Have anesthesia give & hold a deep inspiration to lift the heart off of the spine;
  - 3) Advance primary curve behind the LA. (If no CT, anticipate this is the situation.)
- \*If Valsalva fails, put primary curve to desired side behind trachea & advance ~1 cm q 5 sec.

**Contraindications** are the same as for inserting an OG tube or TEE.

- Upper GI/esophageal disease or abnormalities, including a Hx of:
- Surgery, stricture, varices, hematomas, bleeding, tumors, severe hiatal hernia, etc.

The lung may reduce deflection. Curve doesn't reach spinal border.




### If deflection crosses mid-spine but doesn't reach the spinal border, use Apnea

- 1) Ask Anesthesia to go apneic for 15-20 seconds to deflate the lungs;
- 2) Gently slide stylet ~4" in and out of OG Tube 2-4x to work the curve between spine & lung;
- 3) Position curve behind LA across from Rx area.

### There are 3 EsoSure positions, and 2 ways to change from Right deflection to Left deflection

#### #1- R Deflection Primary Curve




#### #2- L Deflection with Secondary Curve




- a. Advance the Temp Probe into stomach.
- b. Advance OGT/ EsoSure for L deflection with secondary curve.
- c. Then retract Temp Probe.

#### #3- L Deflection with Primary Curve



To use the larger primary curve for L deflection rotate the stylet:

- a. Advance TP.
- b. Retract the EsoSure from the OGT so the curve is above the CS.



- c. Rotate handle 2-3 times, then slowly retract ~4" out of OGT, as rotation usually occurs above the heart. Watch with Fluoro. If it spins 360, rotate in the opposite direction.

- d. After rotation to correct side, release torque & advance curve to desired position.


X When the primary curve is below the CS to the right, do not rotate the EsoSure, only retract it.

X When the primary curve is to the left, do not advance the EsoSure past diaphragm, only rotate.

- e. Then retract TP & position inside stylet's curve.

### Safety Tips

- If resistance is felt while advancing the EsoSure, do not force it. Assess and decide.
- After deflection, visualize the trailing edge of the esophagus using Fluoro with contrast or ICE.
- After deflection, scan LA posterior wall with ICE for indentation on LA caused by displaced anatomy.




E  
q  
u  
i  
p

- **16 or 18 Fr 48" Salem Sump OG Tube (OGT).** Not silicone.
- **20 cc syringe** of room temp water/IV fluid to cool OGT & ...
- **EsoSure** ... **10 cc syringe** for Tube Lube mid-procedure.
- **Tube Lube** is in ea box, for lubricating OGT lumen (Or olive oil).
- **Esophageal Temp Probe (TP).** Use a smooth shaft 9-12 Fr model and avoid Acoustascope models with a balloon over the thermistor/tip.
- **Peds 4-5 mm ET tube or nasal trumpet airway** to be placed orally as an introducer for 9 Fr TP.

P  
r  
e  
p

- 1) **Ask the patient about Contraindications.** Esophageal/upper GI abnormalities.
- 2) **Evaluate chest CT/MRI,** if available, for a narrow LA to spine distance.
- 3) **Endotracheal Intubation** is preferred for EsoSure use. LMA is contraindicated.
- 4) **Lubricate OGT lumen 2 ways:** 1) If EsoSure use is known, suction Tube Lube through OGT from packet before OGT insertion. 2) If EsoSure use is decided mid-procedure, use 10 cc syringe to inject through OGT when use confirmed.
- 5) **Lubricate TP & OGT shaft** with Surgilube. Insert OGT & TP before Heparin. Use a Peds 4-5 mm ET or nasal trumpet (orally) as an introducer for 9 Fr TPs.
- 6) **Insert OG Tube prior to anticoagulation** and remove if EsoSure is not used. Confirm OGT & TP are in the stomach with Fluoro. Insert OG Tube to ~55 cm.
- 7) **Evaluate baseline position of esophagus** to plan deflection work flow.
- 8) **Inject contrast** after transseptal:
  - a) Verify OGT in stomach  $\bar{c}$  Fluoro.
  - b) Retract OGT tip just below CS.
  - c) Inject 10-20 cc contrast slowly.
  - d) Adjust OGT up/down as needed.
  - e) Re-advance OGT into stomach.
 Contrast is not required but is recommended.



I  
n  
s  
e  
r  
t


**1- Bolus Propofol or sedatives** to relax airway & avoid a gag reflex. Wiggle ET tube if unsure. **Inject Tube Lube** in OGT if not already lubricated.

**2- Open & align airway,** in sniffing position by a second person. A two handed head tilt-jaw lift will align oral-esophageal axis for easier insertion.


**3- Cool OGT** to slow stylet's curving, by injecting 20 cc room temp IV fluid/water through OGT connector over 5-10 sec. Hold OGT so it is angled up 45° & straight from mouth.



**4- Advance EsoSure 2-3"** at a time, immediately after injecting cool fluid. Keep OGT straight & watch with Fluoro past throat. \* Advance OGT & EsoSure together if tip stops above CS. Stop when EsoSure tip is past CS or primary curve is at upper heart.




**5- Rotate primary curve** to the desired side. With the primary curve behind the heart, rotate the handle 2-3 times and slowly retract ~4" out of OGT. Rotation usually occurs behind the trachea. Then advance the primary curve back down behind the left atrium.



\* If stylet stops advancing, the most common cause is a loop in OGT. To resolve, retract stylet ~2", fix stylet position and retract OGT over the stylet 2" 3X, then re-advance the stylet. If it won't advance look with Fluoro & remove if no loop.

**6- Initiate a Valsalva** if esophagus does not cross mid-spine: Rotate primary curve to the desired side behind the trachea, then during a Valsalva advance the curve behind the heart. If deflection passes mid-spine but is minimal: Go Apneic for ~20 sec, then slide EsoSure 4" in and out of OGT 2-3 times.

**7- Temp Probe is positioned** inside maximum EsoSure curve across from treatment area.



**8- Assess the esophageal trailing edge** with Fluoro & contrast, or US, while sliding EsoSure in and out to wiggle esophagus. Last, use US to check LA posterior wall for possible indentation.

**EsoSure Removal:**

- Keeping the stylet tip above the diaphragm, hold the end of the OGT and smoothly slide the EsoSure out and in ~4" 2X to slide the esophagus R & L behind the LA. It is hypothesized that side to side esophageal movement may break potential thermal adhesions between tissue layers.
- Lastly, pull EsoSure out of the OGT and suction the stomach & esophagus, especially if contrast was used.